
Pulsed Field Magnetometry Dual Pulse PFM 02

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Hirst

Hirst Magnetic Instruments Ltd

In the continuing search for new Permanent Magnet materials, researchers need to measure magnetic characteristics at higher and higher fields using shorter measurement cycle times.

Pulsed Field Magnetometry not only offers an alternative to existing techniques but also offers major advantages over all established methods.

This booklet provides information on the Pulsed Field Magnetometer for research application developed by Hirst Magnetic Instruments Ltd.

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
Introduction

With the arrival and continuing development of high-coercivity rare-earth permanent magnets, increasingly higher magnetic fields need to be applied to these materials to measure their magnetic characteristics.

The limited fields of electromagnets (2-2.4 Tesla) has led to the widespread use of superconducting magnets in Vibrating Sample Magnetometers (VSM's). These are typically limited to fields of less than 10 Tesla and have the inconvenience of requiring liquid Helium cooling. Their slow performance (measurement cycles of 10 minutes plus) make them unsuitable for many applications and a time consuming operation for others. Where superconducting magnets are required to generate higher fields their costs become prohibitive.

Academic and research institutes developing new materials need to apply high fields of 20 Tesla plus to

determine their full magnetic characteristics including anisotropy data obtained through Single Point Detection (SPD) techniques. This can be achieved utilising Pulsed Field Magnetometry without the need for liquid Helium cooling and with considerably shorter measurement cycle times than available with other techniques.

Pulsed Field Magnetometry

Pulsed Field Magnetometry utilises capacitive discharge to generate high magnetic fields - current laboratory record of 60 Tesla - in conventional resistive solenoids. The duration of the magnetising current is very short, of the order of milliseconds.

The technique therefore employs short measurement cycles with a small time overhead of post processing of data giving overall measurement cycles significantly shorter than available through other techniques.

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
Simulation using finite element analysis of Dynamic magnetisation of NdFeB during magnetisation pulse.

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PFM Theory

Samples of permanent magnetic material to be measured, in cylindrical or similar regular shape, are exposed to a large applied pulsed magnetic field. Inductive pick up coil arrangements are positioned so that they are magnetically coupled with the sample under measurement.

The inductive pick up coil arrangement detects the applied magnetic field (H) and the magnetisation of the sample (M) the signals are digitised by high speed Analogue to Digital Converters (ADC's) and stored in memory.

This recorded raw data is then processed to remove eddy current effects, self demagnetisation factors etc and obtain the four quadrant J/H hysteresis loop measurement, anisotropy data etc.

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Dual Pulse Technique *

A significant innovation offered only by Hirst Magnetic Instruments Ltd, is that of mathematical reduction of eddy current errors.

When a conductive material is subjected to a changing magnetic field, eddy currents will be induced, introducing possible errors in the measurement.

The PFM can be configured to keep these errors small. The Dual Pulse Technique can dramatically reduce any eddy current errors or confirm that they are insignificant.

Applying a magnetic field of a known rate of change to a sample will result in an unknown eddy current and unknown magnetisation of the sample. Application of a magnetic field with a different but known rate of change will result in a different eddy current but the same magnetisation. As the rates of

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change of magnetic fields are known, the relationship between the resultant magnetic fields are known and can therefore be calculated and removed from the magnetisation during the automatic post processing. This techniques is unique to Hirst Magnetic Instruments Ltd.

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
Hirst Magnetic Instrument's Dual Pulse, Pulsed Field Magnetometer PFM02

PFM02


Full loop measurements of sintered and bonded Neodymium Iron Boron and Barium Ferrite.

Hirst Magnetic Instruments' 20 Tesla Dual Pulse, Pulsed Field Magnetometer PFM02 with manual sample loading.

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Instrumentation Unit
Das Meßinstrument

The PFM02 Pulsed Field Magnetometer utilises the Dual Pulse Technique and is designed to generate fields up to 20 Tesla for the measurement of Permanent Magnet materials for research and development applications.

The Capacitive Discharge Unit controls, stores and discharges the electrical charge on the main capacitor bank. When discharged current flows from this unit through the Measurement Unit.

The Measurement Unit contains the main magnetising solenoid, inductive pick up coils, magnet sample and sample holding mechanism together with temperature controlled environment (option). During the discharge signals from the inductive pick up coil arrangement are fed to the Instrumentation Unit.

The Instrumentation Unit integrates and digitises the incoming measurement signals and caches them in a high speed memory. This is then down loaded after the discharge to the host computer for post processing, data storage and record handling.

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Samples to be tested are mounted in the sample holders and are then individually attached to the plunger and inserted into the Measurement Chamber of the PFM02 for measurement (manual system). Using the host PC the researcher initiates pre programmed measurement routines.

The PFM02 then carries out an automatic process, provided that all safety interlocks are correct.

Data is downloaded to the Host PC and the researcher permits the PFM02 to process the raw data by a pre programmed post processing routine to give full magnetic measurement information. Records of measurements are filed on the PC hard disk in a structured form including comprehensive data.

The results are displayed in graphical and tabular form on the Host PC with hard copy being supplied by colour ink jet printer or plotter (option).

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The pre programming of measurements and post processing routines enables a high throughput of samples. All pre programming is carried out by selecting options within a windowing user interface.

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
Various sizes of sample can be positioned in the appropriate sample holder. These are fitted to the sample rod and inserted into the Measurement Chamber. (Manual sample system)

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The PFM02 magnetising solenoid assembly and suspension with Temperature Controlled Environment removed.

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Option - Temperature Controlled Environment

The PFM02 can be supplied with a Temperature Controlled Environment in which samples can be measured at temperatures other than room temperature. This is particularly important when considering Anisotropy measurements. For details contact Sales Office.

Option - Automatic Sample Loading

Enabling multiple samples to be measured automatically to pre programmed requirements. For details contact our Sales Office.

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Option - Colour printer

A suitable ink jet printer/plotter can be supplied to produce colour hard copy. For detailed specification contact our Sales Office.

Option - Additional pickup coils

These are available for smaller sample sizes including flake samples. For details contact Sales Office.

Option - Sample Holders

Sample Holders are supplied with the PFM02. Additional holders can be supplied or easily manufactured locally.

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PFM02 inductive pickup coil assembly and adjustment mechanism.

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

Vacodym 370HR sintered Neodymium Iron Boron. 15 T full loop and second quadrant J and B detail.

PFM02 is a *Dual Pulse Pulsed Field Magnetometer* designed to generate 20 Tesla pulsed fields for the Magnetic Materials Research and Development application specification:

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Hirst Magnetic Instruments Ltd reserve the right to change specifications without notice.

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Hirst Magnetic Instruments Ltd also manufactures wide ranges of magnetic instruments, magnetisers, demagnetisers, precision demagnetisers and special magnetic systems.

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